

Custom Fitting INSTRUCTIONS

"Since 1996"



Give More Care for LESS COST ... The LEEDer way

Fitting Instructions: Multi Use Boot. LMUA, P

Custom fit without tools or heat. Bend it as you see fit.

New

**Softest Adjustable
Interface Available
Flex-E-Core
Technology**

Mold Boot by
HAND as
Necessary to
Custom Fit: Heel
Float, Plantar
Fasciitis,
Night Splint,
Inversion,
Everson

MULTI USE
BOOT /
NIGHT SPLINT

Dorsiflexion
Adjustment
Straps. One
on each side.

Foot Cover
Adjustment
Strap

Walking Sole
Simply wipe clean.

Instep
Strap



**NOW with
Dorsi-Wedge
Dorsiflexion
Pad for
EXTRA Toe
STRETCH
Included
Free!**

**To Adjust: Open all Velcro Closures and Adjust to
Fit Condition. Hand Bend as needed without tools or heat.**

Adult/Pediatric Sizes plus "FIRM"

LMUA_{adult} add NS for Night Splint

LMUP_{pediatric} add NS for Night Splint

Fitting Instructions: Multi Use Boot. LMUA, P

Toll Free Customer Service 866-814-0192

www.leedergroup.com

Toll Free FAX 866-818-0373

leedergroup@leedergroup.com

8508 NW 66th St Miami FL 33166

Custom Fitting INSTRUCTIONS “Since 1996”



Give More Care for LESS COST ... The LEEDer way

Congratulations: you've chosen a product with LEEDer's Flex-E-Core technology: the Multi Use Boot. There is no other Adjustable boot product like this on the market. Custom fit without tools or heat. Bend to fit.

Product Codes: **LMUA**adult, **LMUP**ediatric ADD NS for Night Splint

Make the Multi Use Boot fit perfectly, BEND IT TO FIT. Custom fitting is EASY:

1. If the FOOT section is too long or too short, disconnect the side (dorsiflexion) adjustment straps and adjust the point of bend above or below the heel to shorten or lengthen the splint.

2. Then, RE-BEND the boot so that the FOOT section fits the user, CUSTOM FIT it either longer or shorter.

3. Inversion or Eversion: twist the whole frame to match and correct the problem.

4. Night Splinting: wear with sock; adjust foot position first then “lock” in place with the dorsiflexion straps to counteract plantar flexion contracture; apply pressure in small increments. Cut off straps as desired. OPTION: After fully adjusted, install foam dorsiflexion wedge for extra toe Stretch.

5. For HEEL FLOAT make small bends along the “heel” surface of the boot. Make small bends and small changes to get the perfect fit.

6. For Best Bending Results: USE a “rounded” edge, such as the side rails of a bed or top of the backrest to a chair.

Cut off straps as needed.

Toll Free Customer Service 866-814-0192
www.leedergroup.com

Toll Free FAX 866-818-0373
leedergroup@leedergroup.com

8508 NW 66th St Miami FL 33166